

Product Datasheet

CD81 Antibody (orb1821374)

Description	CD81 Monoclonal Antibody [1.3.3.22]
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC
Immunogen	This antibody was raised by immunizing mice with a B-Cell line derived from a Burkitt lymphoma
Target	CD81
Preservatives	Purified antibody (from supernatant) containing PBS + 0.1% sodium azide
Concentration	1mg/ml
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline use at ICC/IF 1:50 – 1:100, Flow Cyt 1:50 – 1:100
Expiration Date	12 months from date of receipt.