

## Product Datasheet

### S100A9 Antibody (orb1821616)

<b>Description</b>	S100A9 Antibody
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IHC, WB
<b>Immunogen</b>	The original antibody was raised by immunizing Balb/c mice with human blood cells.
<b>Target</b>	S100A9
<b>Preservatives</b>	PBS with 0.02% Proclin 300.
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
<b>Note</b>	For research use only
<b>Application notes</b>	<p>This antibody was used for ELISA on recombinant S100A9 protein (CN104894074). Western blot was performed on S100A9, S100A8, S100A12 and S100A13 using this antibody (CN104894074).</p> <p>Immunohistochemistry was performed on human breast cancer, colon cancer, and hepatocellular carcinoma tissues using this antibody (CN104894074).</p>
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<b>P06702</b>
<b>Expiration Date</b>	12 months from date of receipt.