

Product Datasheet

Biotin Antibody (orb1821627)

Description	Biotin Antibody
Species/Host	Mouse
Reactivity	All
Conjugation	Unconjugated
Tested Applications	ELISA, Epitope Tagging, WB
Immunogen	The original antibody was generated by immunizing BALB/c mice with biotin tagged protein.
Target	Biotin
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	This antibody is recommended for detection of biotin labelled proteins in a sample using ELISA and western blot.
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.