

Product Datasheet

Cytochrome P450 Reductase antibody [Out of stock] (orb182171)

Description

Rabbit polyclonal antibody against

conjugated to Alkaline Phosphatase.

Species/Host

Rabbit

Conjugation

AP

Tested

ICC, IHC, WB

Applications
Immunogen

Rat native full-length Cytochrome P450 Reductase purified from liver tissue

Target

Cytochrome P450 Reductase

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

1 µg/ml of SPC-210 was sufficient for detection of Cytochrome P450 in 10 µg of rat brain lysate by colorimetric immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.

Clonality

Polyclonal

MW

77kDa

Uniprot ID
P00388
NCBI
NP_113764.1
Entrez
29441
Dilution Range

WB (1:1000), ICC/IF (1:60)

Expiration Date

12 months from date of receipt.