

## Product Datasheet

### Cytochrome P450 Reductase antibody [Out of stock] (orb182178)

**Description**

Rabbit polyclonal antibody against

conjugated to RPE.

**Species/Host**

Rabbit

**Conjugation**

RPE

**Tested**

ICC, IHC, WB

**Applications**
**Immunogen**

Rat native full-length Cytochrome P450 Reductase purified from liver tissue

**Target**

Cytochrome P450 Reductase

**Preservatives**

PBS pH 7.4, 50% glycerol, 0.09% sodium azide \*Storage buffer may change when conjugated

**Concentration**

1 mg/ml

**Storage**

Conjugated antibodies should be stored at 4°C

**Note**

For research use only

**Application notes**

1 µg/ml of SPC-210 was sufficient for detection of Cytochrome P450 in 10 µg of rat brain lysate by colorimetric immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.

**Clonality**

Polyclonal

**MW**

77kDa

**Uniprot ID**
**P00388**
**NCBI**
**NP\_113764.1**
**Entrez**
**29441**
**Dilution Range**

WB (1:1000), ICC/IF (1:60)

**Expiration Date**

12 months from date of receipt.