

Product Datasheet

Anti-VMAT2 Antibody (orb1821795)

Description	Anti-VMAT2 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide corresponding to amino acids from an internal sequence of the human VMAT2 protein
Preservatives	PBS + 0.02% NaN ₃
Form/Appearance	Antigen Affinity Purified
Storage	Recommended that the undiluted antibody be aliquoted into smaller working volumes (10-30 uL/vial depending on usage) upon arrival and stored long term at -20° C or -80° C, while keeping a working aliquot stored at 4° C for short term. Avoid freeze/thaw cycles. Stable for at least 1 year.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q05940
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.