

Product Datasheet

Anti-Antennapedia Homeodomain Antibody (orb1821879)

Description

Anti-Antennapedia Homeodomain

Species/Host	Rabbit
Reactivity	All
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide corresponding to amino acids from the Antennapedia homeodomain sequence which is 100 % conserved across all known species.
Preservatives	PBS + 0.02% NaN ₃
Form/Appearance	Antigen Affinity Purified
Storage	Recommended that the undiluted antibody be aliquoted into smaller working volumes (10-30 uL/vial depending on usage) upon arrival and stored long term at -20° C or -80° C, while keeping a working aliquot stored at 4° C for short term. Avoid freeze/thaw cycles. Stable for at least 1 year.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P02833
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.