

Product Datasheet

Anti-N-WASP Antibody (orb1821943)

Description	Anti-N-WASP Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-WASP synthetic peptide (coupled to KLH) corresponding to amino acid residues in the N-terminal region of human N-WASP. This N-WASP peptide sequence is 100% homologous to rat and mouse N-WASP, and has low homology to the corresponding region in the human WASP.
Preservatives	PBS + 1 mg/ml rAlbumin, 0.05% NaN ₃ and 50% glycerol
Form/Appearance	Antigen Affinity Purified
Storage	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	65
Uniprot ID	O00401
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.