

Product Datasheet

Anti-WASP / N-WASP Antibody (orb1821944)

Description	Anti-WASP / N-WASP Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	A synthetic peptide (coupled to KLH) corresponding to a region near the C-terminus of human WASP. This sequence is highly homologous with a similar region in human N-WASP.
Preservatives	PBS + 1 mg/ml rAlbumin, 0.05% NaN ₃ and 50% glycerol
Form/Appearance	Antigen Affinity Purified
Storage	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	63
Uniprot ID	O00401
Dilution Range	WB: 1:1000, ICC: 1:100
Expiration Date	12 months from date of receipt.