

Product Datasheet

HO-1 antibody (orb182195)

Description	Rabbit polyclonal antibody against HO-1, conjugated to PerCP.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PerCP
Tested Applications	IHC, IP, WB
Immunogen	Rat native full-length HO-1 purified from liver tissue
Target	HO-1
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	5 µg/ml of SPC-211 was sufficient for detection of HO-1 in 10 µg of rat brain lysate by colorimetric immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	33kDa
Uniprot ID	P06762
NCBI	NP_036712.1
Entrez	24451
Dilution Range	WB (1:400), IHC (1:1000), IP (1:100)
Expiration Date	12 months from date of receipt.