

Product Datasheet

Anti-Trop-2/TACSTD2 (Extracellular region) Antibody (orb1821973)

Description	Anti-Trop-2/TACSTD2 (Extracellular region) Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IP, WB
Immunogen	Clone M005 was generated from a proprietary antigen related to the extracellular region of human Trop-2 from the A431 epidermoid carcinoma cell line.
Preservatives	PBS + 1 mg/ml rAlbumin, 0.05% NaN3 and 50% glycerol
Form/Appearance	Protein G Purified
Storage	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	50-65
Uniprot ID	P09758
Dilution Range	WB: 1:1000, ICC: 1:100, IP: 1:100
Expiration Date	12 months from date of receipt.