

Product Datasheet

Gelsolin Gsn Antibody (Monoclonal, GS-2C4) (orb18221)

Description

Gelsolin Gsn Antibody (Monoclonal, GS-2C4)

Species/Host

Mouse

Reactivity

Human, Rabbit

Conjugation

Unconjugated

Tested

IHC, WB

Applications

Immunogen

Human plasma gelsolin.

Form/Appearance

Lyophilized

Concentration

Adding 1 ml of PBS buffer will yield a concentration of 100 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 1ml of PBS buffer will yield a concentration of 100ug/ml.

Isotype

Mouse IgG1

Clonality

Monoclonal

MW

86068 MW

Uniprot ID


[Q68FP1](#)

Dilution Range


Immunohistochemistry (Paraffin-embedded Section), 1-2µg/ml, Human, rabbit, By Heat
Immunohistochemistry (Frozen Section), 1-2µg/ml, Human, rabbit, -Western blot, 0.5-1µg/ml, Human, rabbit

Expiration Date


12 months from date of receipt.



IHC(P)
analysis of
Human
Lung
Cancer Tis...



IHC(P)
analysis of
Human
Intestinal
Canc...



WB analysis
of Gelsolin
using anti-
Gelso...