

Product Datasheet

Anti-Laminin β 2/ γ 1 Subunits Antibody (orb1822130)

Description	Anti-Laminin β 2/ γ 1 Subunits Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Clone (M046) was generated from a recombinant human laminin 521 protein that includes a heterotrimer of three subunits laminin α 5, laminin β 2, and laminin γ 1 subunits.
Preservatives	PBS + 1 mg/ml rAlbumin, 0.05% NaN ₃ and 50% glycerol
Form/Appearance	Protein A Purified
Storage	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	760
Uniprot ID	P11047
Dilution Range	WB: 1:1000, ICC: 1:200
Expiration Date	12 months from date of receipt.