

Product Datasheet

Anti-EGFR (Ser-967), Phosphospecific Antibody (orb1822191)

Description

Anti-EGFR (Ser-967), Phosphospecific

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Phospho-EGFR (Ser-967) synthetic peptide (coupled to KLH) corresponding to amino acid residues surrounding serine 967 of human EGFR (ErbB-1). This human EGFR sequence has high homology to rat and mouse EGFR, and low homology to other ErbB family members.
Preservatives	PBS + 1 mg/ml rAlbumin, 0.05% NaN ₃ and 50% glycerol
Form/Appearance	Antigen Affinity Purified
Storage	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	180
Uniprot ID	P00533
Dilution Range	WB: 1:1000, ICC: 1:200
Expiration Date	12 months from date of receipt.