

Product Datasheet

Anti-eEF2K (Ser-445), Phosphospecific Antibody (orb1822195)

Description

Anti-eEF2K (Ser-445),

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat, Xenopus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Phospho-eEF2K (Ser-445) synthetic peptide (coupled to carrier protein) corresponds to amino acids surrounding Ser-445 in human eEF2K. This sequence is well conserved in rat, mouse, and fish eEF2K, and has low homology to other proteins.
Preservatives	PBS + 1 mg/ml rAlbumin, 0.05% NaN3 and 50% glycerol
Form/Appearance	Antigen Affinity Purified
Storage	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	105
Uniprot ID	000418
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.