

Product Datasheet

Anti-CD63 (Extracellular region) Antibody (orb1822260)

Description	Anti-CD63 (Extracellular region) Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Clone (M040) was generated from a proprietary antigen related to the native human CD63 expressed in MeWo melanoma cell line.
Preservatives	PBS + 1 mg/ml rAlbumin, 0.05% NaN ₃ and 50% glycerol
Form/Appearance	Protein G Purified
Storage	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	30-60
Uniprot ID	P08962
Dilution Range	WB: 1:100, ICC: 1:500
Expiration Date	12 months from date of receipt.