

Product Datasheet

Anti-E-Cadherin (Cytoplasmic) Antibody (orb1822306)

Description

Anti-E-Cadherin (Cytoplasmic)

Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, IP, WB
Immunogen	Clone (M168) was generated from a mouse recombinant E-Cadherin protein containing amino acids in the cytoplasmic region. This sequence is highly conserved in human and rat E-cadherin.
Preservatives	PBS + 1 mg/ml rAlbumin, 0.05% NaN ₃ and 50% glycerol
Form/Appearance	Protein A Purified
Storage	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	120
Uniprot ID	P09803
Dilution Range	WB: 1:1000, IHC: 1:50, ICC: 1:250, IP: 1:100
Expiration Date	12 months from date of receipt.