

Product Datasheet

PDI Antibody (orb1822442)

Description	Rabbit Anti-Rat PDI Polyclonal
Species/Host	Rabbit
Reactivity	Bovine, Canine, Frog, Guinea pig, Hamster, Human, Invertebrate, Mouse, Mussel, Porcine, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	ICC, IHC, IP, WB
Immunogen	Rat PDI synthetic peptide (409-509 aa C-terminal) conjugated to KLH
Target	PDI
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	-20°C
Note	For research use only
Application notes	A 1:1000 dilution of SPC-114 was sufficient for detection of PDI in 20 µg of HeLa cell lysate by ECL immunoblot analysis.
Clonality	Polyclonal
MW	58kDa
Uniprot ID	P04785
NCBI	NP_001099245.2
Entrez	287164
Dilution Range	WB (1:1000), ICC/IF (1:100); optimal dilutions for assays should be determined by the user.
Expiration Date	12 months from date of receipt.