

Product Datasheet

ACE2 Antibody / Angiotensin-Converting Enzyme 2 (orb1822584)

Description

Carboxypeptidase which converts

peptide of unknown function, and...

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC-P

Applications
Immunogen

A recombinant partial protein sequence (within amino acids 705-805) from the human protein was used as the immunogen for the Angiotensin Converting Enzyme 2 antibody.

Storage

Aliquot the Angiotensin Converting Enzyme 2 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

Note

For research use only

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Rabbit IgG, kappa

Clonality

Recombinant

Uniprot ID

Q9BYF1

Hazard
Information

This Angiotensin Converting Enzyme 2 antibody is available for research use only.

Dilution Range

ELISA (Order BSA-free format for coating), Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.