

Product Datasheet

S100G Antibody / Calbindin D9K (orb1822599)

Description	This gene encodes calbindin D9K, a vitamin D-dependent calcium-binding protein. This cytosolic prote...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein sequence (within amino acids 1-79) from the human protein was used as the immunogen for the S100G antibody.
Storage	Aliquot the S100G antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2, kappa
Clonality	Monoclonal
Uniprot ID	P29377
Hazard Information	This S100G antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.