

Product Datasheet

Mouse IgG2a Isotype Control Antibody (orb1822665)

Description	To serve as a true negative control for various procedures, this rabbit monoclonal IgG should be use...
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FACS, IF, IHC-P, WB
Immunogen	Recombinant full-length human IG1R protein was used as the immunogen for the Mouse IgG2a Isotype Control antibody.
Storage	Aliquot the Mouse IgG2a Isotype Control antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2a
Clonality	Monoclonal
Hazard Information	This Mouse IgG2a Isotype Control antibody is available for research use only.
Dilution Range	ELISA: assay dependent,Flow Cytometry: assay dependent,Immunofluorescence: assay dependent,Western blot: assay dependent,Immunohistochemistry (FFPE): assay dependent
Expiration Date	12 months from date of receipt.