

Product Datasheet

SNAPC4 Antibody (orb1822678)

Description	TATA-box binding protein (TBP) interactions with TBP-associated factors (TAFs) are required for the ...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FACS
Immunogen	Recombinant full-length human protein was used as the immunogen for the SNAPC4 antibody.
Storage	Aliquot the SNAPC4 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2b, kappa
Clonality	Monoclonal
Uniprot ID	Q5SXM2
Hazard Information	This SNAPC4 antibody is available for research use only.
Dilution Range	Flow cytometry: 1-2ug/million cells
Expiration Date	12 months from date of receipt.