

Product Datasheet

TrkB Antibody / NTRK2 (orb1822725)

Description	The Trk proto-oncogene encodes a tyrosine protein kinase, Trk A, also designated Trk gp140, that ser...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein sequence (within amino acids 250-450) from the human protein was used as the immunogen for the Tropomyosin-related kinase B antibody.
Storage	Aliquot the Tropomyosin-related kinase B antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2c, lambda
Clonality	Monoclonal
Uniprot ID	Q16620
Hazard Information	This Tropomyosin-related kinase B antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.