

Product Datasheet

PD-L1 Antibody / B7-H1 / CD274 (orb1822757)

Description	PD-L1 plays a critical role in induction and maintenance of immune tolerance to self. [UniProt]
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein sequence (within amino acids 190-290) from the human protein was used as the immunogen for the Programmed cell death 1 ligand 1 antibody.
Storage	Aliquot the Programmed cell death 1 ligand 1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	Q9NZQ7
Hazard Information	This Programmed cell death 1 ligand 1 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.