

Product Datasheet

NECAB1 Antibody / N-terminal EF-hand calcium binding protein 1 (orb1822858)

Description	Members of the EF-CBP (N-terminal EF-hand calcium binding protein)/NECAB (neuronal calcium-binding p...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant human NECAB1 protein was used as the immunogen for the NECAB1 antibody.
Storage	Aliquot the NECAB1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG
Clonality	Monoclonal
Uniprot ID	Q8N987
Hazard Information	This NECAB1 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2 µg/ml for 30 min at RT
Expiration Date	12 months from date of receipt.