

Product Datasheet

NECAB1 Antibody / N-terminal EF-hand calcium binding protein 1 (orb1822861)

Description

Members of the EF-CBP (N-terminal

protein)/NECAB (neuronal calcium-binding p...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications
Immunogen

Recombinant human NECAB1 protein was used as the immunogen for the N-terminal EF-hand calcium binding protein 1 antibody.

Storage

Aliquot the N-terminal EF-hand calcium binding protein 1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

Note

For research use only

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG

Clonality

Monoclonal

Uniprot ID
Q8N987
Hazard
Information

This N-terminal EF-hand calcium binding protein 1 antibody is available for research use only.

Dilution Range

Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.