

Product Datasheet

Insulin Ins2 Antibody (Monoclonal, K36AC10) (orb18229)

Description

Insulin Ins2 Antibody (Monoclonal, K36AC10)



Species/Host

Mouse

Reactivity

Human, Rabbit, Rat


Conjugation

Unconjugated

Tested

IHC

Applications



Immunohistochemical analysis of paraffin...

Immunogen

Human insulin.

Form/Appearance

Lyophilized

Concentration

Adding 1 ml of PBS buffer will yield a concentration of 100 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 1ml of PBS buffer will yield a concentration of 100ug/ml.

Isotype

Mouse IgG1

Clonality

Monoclonal

MW

12339 MW

Uniprot ID

P01323

Dilution Range

Immunohistochemistry (Paraffin-embedded Section), 0.4-1µg/ml, Human, rabbit, rat, By Heat

Expiration Date

12 months from date of receipt.