

Product Datasheet

WT1-associated protein Antibody / WTAP (orb1823037)

Description	Wilms tumor (WT) is an embryonal malignancy of the kidney that affects 1 in 10,000 infants and is ob...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FACS, IHC-P
Immunogen	Recombinant full-length human protein was used as the immunogen for the WT1-associated protein antibody.
Storage	Aliquot the WT1-associated protein antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2b, kappa
Clonality	Monoclonal
Uniprot ID	Q15007
Hazard Information	This WT1-associated protein antibody is available for research use only.
Dilution Range	Flow cytometry: 1-2ug/million cells, Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.