

Product Datasheet

Mineralocorticoid Receptor Antibody / MR / NR3C2 (orb1823379)

Description	Mineralocorticoid hormones are primarily found in epithelial tissues where they function as regulato...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FACS
Immunogen	A recombinant partial protein sequence (within amino acids 601-673) from the human protein was used as the immunogen for the NR3C2 antibody.
Storage	Aliquot the NR3C2 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1
Clonality	Monoclonal
Uniprot ID	P08235
Hazard Information	This NR3C2 antibody is available for research use only.
Dilution Range	Flow cytometry: 1-2ug/million cells
Expiration Date	12 months from date of receipt.