

Product Datasheet

IL-2 Antibody (orb1823411)

Description	Interleukin-2 (IL-2) is a T cell stimulatory cytokine that induces T cell proliferation and NK cell ...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FACS
Immunogen	A recombinant fragment of human IL2 protein was used as the immunogen for the Interleukin-2 antibody.
Storage	Aliquot the Interleukin-2 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Mouse IgG1
Clonality	Monoclonal
Uniprot ID	P60568
Hazard Information	This Interleukin-2 antibody is available for research use only.
Dilution Range	ELISA (Order BSA-free format for coating),Flow cytometry: 1-2ug/million cells
Expiration Date	12 months from date of receipt.