


Product Datasheet

Mouse Amphiregulin(AR) ELISA Kit (orb182343)

Description	Mouse Amphiregulin(AR) ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: V100-A248
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	27549 MW
Uniprot ID	P31955
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Mouse
Amphiregulin/AR
PicoKine ELISA
Kit...