


Product Datasheet

Mouse VE-Cadherin Cdh5 ELISA Kit (orb182346)

Description	Mouse VE-Cadherin Cdh5 ELISA Kit
Reactivity	Mouse
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: D46-Q599
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	87903 MW
Uniprot ID	P55284
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:100, concentration in serum and plasma is 100-200 ng/ml.
Expiration Date	Please enquire.



Mouse VE-Cadherin/CD144 PicoKine ELISA Kit