


Product Datasheet

Human CD93/C1qR ELISA Kit (orb182347)

Description	Human CD93/C1qR ELISA Kit
Reactivity	Human
Range	78 pg/ml - 5,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A24-K580
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA), saliva and urine.
MW	68560 MW
Uniprot ID	Q9NPY3
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:100, concentration in serum and plasma is 72-267 ng/ml.
Expiration Date	Please enquire.



Human
CD93/C1qR
PicoKine
ELISA Kit
stand...