

Product Datasheet

CD10 Antibody / Neprilysin (orb1823568)

Description Recognizes a 100kDa glycoprotein, identified as CD10, also known as Common Acute Lymphocytic Leukemi...

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications IHC-P, WB

Immunogen A recombinant human CD10 protein fragment (within amino acids 1-200) was used as the immunogen for the CD10 antibody.

Storage Aliquot the CD10 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

Note For research use only

Formula 0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype Mouse IgG1, kappa

Clonality Monoclonal

Uniprot ID **P08473**

Hazard Information This CD10 antibody is available for research use only.

Dilution Range Western blot: 1-2ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

Expiration Date 12 months from date of receipt.