

Product Datasheet

I-309/CCL1 IGF1R Antibody (orb182360)

Description

I-309/CCL1 IGF1R Antibody

97KD — 1 2 3

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

WB analysis
of I-309
using anti-I-
309 an...

Immunogen

E.coli-derived human I-309 recombinant protein (Position: K24-K96). Human I-309 shares 38% amino acid (aa) sequence identity with mouse I-309.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

WB: The detection limit for I-309 is approximately 0.25ng/lane under reducing conditions. Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

154793 MW

Uniprot ID

P08069

Dilution Range

Western blot, 0.1-0.5µg/ml, Human, -ELISA , 0.1-0.5µg/ml, Human

Expiration Date

12 months from date of receipt.