

Product Datasheet

PDCD1 Antibody / PD-1 / PD1 (orb1823700)

Description	PDCD-1 (programmed cell death-1 protein), also designated CD279, is a type I transmembrane receptor ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	A recombinant partial protein sequence (within amino acids 1-200) from the human protein was used as the immunogen for the PD1 antibody.
Storage	Aliquot the PD1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	Q15116
Hazard Information	This PD1 antibody is available for research use only.
Dilution Range	ELISA (Order BSA-free format for coating),Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.