

Product Datasheet

SMARCAL1 Antibody (orb1823781)

Description	SmrcAL1 (SWI/SNF-related matrix-associated Actin-dependent regulator of chromatin subfamily A-like ...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant full-length human SMARCAL1 protein was used as the immunogen for the SMARCAL1 antibody.
Storage	Aliquot the SMARCAL1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2c, kappa
Clonality	Monoclonal
Uniprot ID	Q9NZC9
Hazard Information	This SMARCAL1 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.