


Product Datasheet


Mucin Gastric MUC5AC Antibody (Monoclonal, 45M1) (orb18238)

Description Mucin Gastric MUC5AC Antibody (Monoclonal,


Species/Host	Mouse
Reactivity	Human, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, WB
Immunogen	Mucin from human ovarian cyst fluid.
Form/Appearance	Lyophilized
Concentration	Adding 1 ml of PBS buffer will yield a concentration of 100 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 1ml of PBS buffer will yield a concentration of 100ug/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	585570 MW
Uniprot ID	P98088
Dilution Range	Immunohistochemistry (Paraffin-embedded Section), 2-4µg/ml, Human, rabbit, rat, By Heat Immunocytochemistry , 1µg/ml, Human, rabbit, rat, -Western blot, 1-2µg/ml, Human, rabbit, rat
Expiration Date	12 months from date of receipt.



IHC
analysis of
Mucin
Gastric
using anti...



IHC
analysis of
Mucin
Gastric
using anti...



IHC
analysis of
Mucin
Gastric
using anti...