

Product Datasheet

PAX8 Antibody (orb1823867)

Description	Recognizes a protein of 62kDa, identified as PAX8. It is a member of the paired box (PAX) family of ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein sequence (within amino acids 150-300) from the human protein was used as the immunogen for the PAX8 antibody.
Storage	Aliquot the PAX8 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	Q06710
Hazard Information	This PAX8 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.