

Product Datasheet

IL18R1 Antibody / Interleukin 18 receptor alpha (orb1823892)

Description

Interleukin-18 (IL-18) has been

IFN-Gamma production and enhan...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications
Immunogen

A recombinant partial protein sequence (within amino acids 1-200) from the human protein was used as the immunogen for the Interleukin 18 receptor alpha antibody.

Storage

Aliquot the Interleukin 18 receptor alpha antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

Note

For research use only

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Uniprot ID

[Q13478](#)

Hazard
Information

This Interleukin 18 receptor alpha antibody is available for research use only.

Dilution Range

Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.