


Product Datasheet

Rab5c antibody (orb182401)

Description	Goat polyclonal antibody to Rab5c. Rab5c belongs to the small GTPase superfamily, Rab family. Rab5c ...
Species/Host	Goat
Reactivity	Bovine, Canine, Donkey, Equine, Feline, Gallus, Goat, Guinea pig, Hamster, Human, Mouse, Other, Porcine, Primate, Rabbit, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	IHC-Fr, IHC-P, WB
Immunogen	Purified recombinant peptide derived from within residues 118 aa to the C-terminus of mouse Rab5c produced in E. coli.
Target	Rab5c
Preservatives	PBS, 20% glycerol and 0.05% sodium azide
Form/Appearance	Polyclonal antibody supplied as a 200 µl (3 mg/ml) aliquot in PBS, 20% glycerol and 0.05% sodium azide. This antibody is epitope-affinity purified from goat antiserum.
Concentration	3 mg/ml
Storage	Store at -20°C for long-term storage. Store at 2-8°C for up to one month.
Note	For research use only
Application notes	The antibody solution should be gently mixed before use.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB:1:250-1:2,000, IHC-P:1:200-1:750, IHC-F:1:200-1:750
Expiration Date	12 months from date of receipt.



Western blot analysis of At-T20 cell lin...