

Product Datasheet

Syndecan-1 Antibody / CD138 (orb1824021)

Description	CD138, also designated syndecan-1, is a member of the syndecan family of four transmembrane spanning...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein sequence (within amino acids 1-250) from the human protein was used as the immunogen for the SCD1 antibody.
Storage	Aliquot the SCD1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Uniprot ID	P18827
Hazard Information	This SCD1 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.