

Product Datasheet

MIF Antibody / Macrophage migration inhibitory factor (orb1824114)

Description	Macrophage migration inhibitory factor, known as MIF or glycosylation inhibiting factor, is a secret...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	A recombinant fragment of human protein was used as the immunogen for the MIF antibody.
Storage	Aliquot the MIF antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG
Clonality	Monoclonal
Uniprot ID	P14174
Hazard Information	This MIF antibody is available for research use only.
Dilution Range	Western blot: 1-2ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.