


Product Datasheet

ADH1 antibody (orb182412)

Description	ADH1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human ADH1A (201-300/375 aa)
Target	ADH1A + ADH1B + ADH1G
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	40 kDa
Uniprot ID	P07327
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
MCF7 cell
US...