

Product Datasheet

Chromogranin B Antibody / Secretogranin I (orb1824128)

Description	Chromogranins (secretogranins) are acidic glycoproteins that localize within secretory granules of e...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein sequence (within amino acids 250-450) from the human protein was used as the immunogen for the Chromogranin B antibody.
Storage	Aliquot the Chromogranin B antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Mouse IgG
Clonality	Monoclonal
Uniprot ID	P05060
Hazard Information	This Chromogranin B antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.