

Product Datasheet

Chromogranin B Antibody / Secretogranin I (orb1824132)

Description	Chromogranins (secretogranins) are acidic glycoproteins that localize within secretory granules of e...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	A recombinant partial protein sequence (within amino acids 250-450) from the human protein was used as the immunogen for the Secretogranin I antibody.
Storage	Aliquot the Secretogranin I antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG, kappa
Clonality	Monoclonal
Uniprot ID	P05060
Hazard Information	This CHGB antibody is available for research use only.
Dilution Range	ELISA (Order BSA-free format for coating)
Expiration Date	12 months from date of receipt.