



Product Datasheet

TAP1 antibody (orb182416)


Description	TAP1 antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human TAP1 (501-600/808 aa)
Target	TAP1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	87 kDa
Uniprot ID	Q03518
Expiration Date	12 months from date of receipt.



IHC-P staining of mouse stomach tissue u...



Immunohistochemical staining of rat ovar...



Western blot analysis of mouse spleen ly...