

Product Datasheet

Chicken IgY antibody (orb182471)

Description	Chicken IgY antibody
Species/Host	Rabbit
Reactivity	Gallus
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	Native Chicken IgY
Preservatives	0.01M PBS (pH7.4)
Form/Appearance	Lyophilized powder
Storage	Store at -20°C stable for 2 years (lyophilized powder). Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,IF=1:100-500,ELISA=1:500-1000
Expiration Date	12 months from date of receipt.