

Product Datasheet

SBSN Antibody / Suprabasin (orb1824720)

Description	Suprabasin (SBSN) is a 247 amino acid secreted protein that is upregulated in differentiating kerati...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant fragment from the human protein was used as the immunogen for the SBSN antibody.
Storage	Aliquot the SBSN antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Uniprot ID	Q6UWP8
Hazard Information	This SBSN antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.