

Product Datasheet

AMH Antibody / Anti-Muellerian Hormone (orb1824812)

Description	The transforming growth factor beta (TGF beta) superfamily is composed of numerous growth and differ...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant fragment of the C-terminus of the human protein (within amino acids 460-560) was used as the immunogen for the Anti Mullerian Hormone antibody.
Storage	Aliquot the Anti Mullerian Hormone antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Uniprot ID	P03971
Hazard Information	This Anti Mullerian Hormone antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.